



CERTIFICATE OF ANALYSIS

Product ANTI-DIGOXIGENIN/DIGOXIN
AFFINITY-PURIFIED

Catalog No. MB-7000

Amount 1 mg

Lot No. Z0214

Produced in goat

Concentration 1 mg/ml

Solution 10 mM phosphate, pH 7.8, 0.15 M NaCl, 0.08% sodium azide.
Sucrose has been added for stabilization.

Storage Refrigerate

Remarks

This product is recommended for use in transfer blot applications (Northern, Southern, Western, dot blots, etc.) and *in situ* hybridization.

The recommended concentration range for use is 0.5-5 µg/ml. This reagent can be detected using a number of options including biotinylated horse anti-goat IgG (Cat. No. BA-9500) or biotinylated rabbit anti-goat IgG (Cat. No. BA-5000) and the VECTASTAIN[®] Elite[®] ABC Kit (Cat. No. PK-6100) or VECTASTAIN[®]-AP ABC Kit (Cat. No. AK-5000) detection system.

This reagent can also be detected using ImmPRESS[™] Polymer Reagent Anti-Goat Ig (MP-7405) or VECTASTAIN[®] Elite ABC Kit, Goat IgG (Cat. No. PK-6105).