



CERTIFICATE OF ANALYSIS

Product **2.5% NORMAL GOAT SERUM**

Catalog No. S-1012

Amount 50 ml

Lot No. ZC1208

Other Remarks:

This serum was collected from normal, healthy animals. The serum was heat-treated and centrifuged to remove precipitates and filtered through a 0.45 μ m filter. This preparation is 2.5% NGS (v/v) diluted in buffer and contains 0.08% sodium azide. This is the same serum preparation that is included in the ImmPRESS™ Detection Systems and the R.T.U. VECTASTAIN® ABC Kits. This reagent can be used for blocking non-specific binding or as an antibody diluent.