



---

## CERTIFICATE OF ANALYSIS

Product     **2.5% NORMAL HORSE SERUM**

Catalog     S-2012

Amount     50 ml

Lot No.     ZC0620

Other Remarks:

This serum was collected from normal, healthy animals. The serum was heat-treated and centrifuged to remove precipitates and filtered through a 0.45  $\mu\text{m}$  filter. This preparation is 2.5% NHS (v/v) diluted in buffer and contains 0.08% sodium azide. This is the same serum preparation that is included in the ImmPRESS™ Detection Systems and the R.T.U. VECTASTAIN® ABC Kits. This reagent can be used for blocking non-specific binding or as an antibody diluent.