



CERTIFICATE OF ANALYSIS

Product AGAROSE-STREPTAVIDIN

Catalog SA-5010

Amount 2 ml of settled gel (1:1, v/v slurry)

Lot No. W1004

Suspension Solution 10 mM HEPES, pH 8.0, 0.15 M NaCl, 0.08% sodium azide.

Protein Concentration 1 mg/ml of settled gel

Binding Capacity ≥1 mole of biotinylated-horseradish peroxidase/mole
of streptavidin

Storage Conditions Refrigerate - DO NOT FREEZE

Other Remarks

Wash gel thoroughly with buffer before use.