



---

## CERTIFICATE OF ANALYSIS

Product DyLight<sup>®</sup> 549 STREPTAVIDIN

Catalog No. SA-5549

Amount 1 mg active conjugate

Lot No. ZC0830

Concentration 1.0 mg/ml active conjugate

Solution 10 mM HEPES, 0.15 M NaCl, pH 7.5, 0.08% sodium azide.

Maximum Excitation @ 556 nm

Maximum Emission @ 571 nm

F/P (molar) 2.8

Storage Conditions Refrigerate in the dark.

Recommended concentration range for use 5-30 µg/ml

Avoid using RPMI 1640 or other biotin-containing solutions as diluents. Serum also can contain biotin and should not be added to diluents.

*DyLight<sup>®</sup> is a trademark of Thermo Fisher Scientific Inc. and its subsidiaries.*