



CERTIFICATE OF ANALYSIS

Product DyLight[®] 649 STREPTAVIDIN

Catalog No. SA-5649

Amount 1 mg active conjugate

Lot No. ZB0424

Concentration 1.0 mg/ml active conjugate

Solution 10 mM HEPES, 0.15 M NaCl, pH 7.5, 0.08% sodium azide.

Maximum Excitation @ 655 nm

Maximum Emission @ 670 nm

F/P (molar) 3.7

Storage Conditions Refrigerate in the dark.

Recommended concentration range for use 5-30 µg/ml

Avoid using RPMI 1640 or other biotin-containing solutions as diluents. Serum also can contain biotin and should not be added to diluents.

DyLight[®] is a trademark of Thermo Fisher Scientific Inc. and its subsidiaries.