



---

## CERTIFICATE OF ANALYSIS

Product **AGAROSE PROTEIN A**

Catalog No. SP-0050

Amount 5 ml of settled gel (1:1, v/v slurry)

Lot No. Z1220

Suspension Solution 10 mM phosphate, pH 7.5, 0.15 M NaCl, 0.08% sodium azide

Protein Concentration 2.0 mg/ml of settled gel

Binding Capacity >10 mg human IgG/ml gel

Storage Conditions Refrigerate - DO NOT FREEZE

### Remarks

Elution can generally be accomplished with several column volumes of 200 mM acetic acid. After use, wash the gel with several column volumes of PBS, then resuspend gel in PBS containing 0.08% sodium azide for storage.