

CERTIFICATE OF ANALYSIS

Product BIOTINYLATED DEXTRAN AMINE (10Kd) FLUORESCEIN BDA-10KF

Catalog No. SP-1130

Amount <u>10 mg</u>

Lot No. ZB2813

Free Fluorochrome _none

Approximate size <u>10,000 MW</u>

Storage Conditions Refrigerate at 4 °C or, after reconstitution, store frozen.

Remarks:

Maximum Excitation @ 495-500 nm Maximum Emission @ 514-521 nm

Reconstitute by addition of 0.4 ml of glass distilled water. The reconstituted product will have the following composition: 15 mM NaCl, pH 8.0. No preservative has been added.