



CERTIFICATE OF ANALYSIS

Product CASEIN SOLUTION (10X)

Catalog No. SP-5020

Amount 250 ml

Lot No. ZB0922

Storage 2-8 °C

Usage Information:

This product is designed to be used as a general blocking agent/diluent for nucleic acid and protein blotting applications. Prior to use, 10x Casein Solution should be diluted to 1x with distilled or deionized water. For prolonged storage of diluted Casein Solution, sodium azide should be added to give a final concentration of 0.08% (w/v) and the diluted solution refrigerated. 1x Casein Solution is the recommended diluent for alkaline phosphatase-conjugated antibodies, streptavidin, and lectins used in the FastTag[®] nucleic acid labeling system.