



## SF2600

The Stat Fax 2600 **Microplate Washer** is fully automatic. Washes flat, round, V-bottom plates and strips. No user adjustments required. Automatic calibration, alignment, and last row detection. Large, non-volatile memory stores at least 50 user wash protocols, precision operation with constant monitoring of vacuum and pressure. Pump cycles only when necessary for quiet operation. Programmable automatic rinse cycle. Ready to run: bottles, aerosol shield, 8-way manifold included. 12 and 16 position wash heads available as accessories. NRTL and CE certified.